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INSTRUCTIONS FOR USE

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TECHNICAL NAME: Box / case for surgical instruments

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TRADE NAME: Cases for Implants and Surgical Instruments AL.

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WARNING: Read carefully all of the instructions before use. Comply with all warnings and precautions mentioned in these instructions. The failure to comply with these recommendations may lead to complications.

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Note: The Cases for Packaging Implants and Surgical Instruments AL are manufactured for packaging implants and / or instrumental used in orthopedic surgeries and bone correction, and for each type of surgery and for each set of products there is a specific case.

Purpose: The purpose of Cases for packaging of surgical implants and instruments is to pack and easy transport and sterilization of instruments and materials - health products - identifying them as a way to facilitate the surgical instruments and sterilization by autoclave.

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2. Instrument Raw Materials


The following materials are employed on manufacturing the cases: Aluminum (AA 6351)

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


These Instructions also regard the following products:

CÓDIGO	DESCRIÇÃO	FUNÇÃO	ESPECIFICAÇÃO TÉCNICA	ILUSTRAÇÃO
PA.02.03.1595	Estojo trauma 1,5 nº1AL Estuche trauma 1.5 nº1 AL	Acomodar implantes e instrumentos	Alumínio (AA 6351)	
PA.02.03.1596	Estojo trauma 1,5 nº2 AL Estuche trauma 1.5 nº2 AL			
PA.02.03.1597	Estojo trauma 1,5 nº3 AL Estuche trauma 1.5 nº3 AL			
PA.02.03.1598	Estojo trauma 2,0 nº1 AL Estuche trauma 2.0 nº1 AL			
PA.02.03.1599	Estojo trauma 2,0 nº2 AL Estuche trauma 2.0 nº2 AL			
PA.02.03.1600	Estojo trauma 2,0 nº3 AL Estuche trauma 2.0 nº3 AL			
PA.02.03.1601	Estojo trauma 2,4 nº1 AL Estuche trauma 2.4 nº1 AL			
PA.02.03.1602	Estojo trauma 2,4 nº2 AL Estuche trauma 2.4 nº2 AL			

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4. Operation principle

Accommodate in a systematic and organized way the implants and instruments in order to maintain the physical integrity of the instruments, facilitating handling both during surgery, as in transportation procedures, storage and sterilization.

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5. Directions

Each model case for packaging surgical implants and instruments has a specific form of packaging of implants and / or instruments that are provided by the IDs each case, facilitating the recognition of the instruments and / or implants. The kits do not come accompanied by instruments and implants for surgeries, these are purchased separately, then put in their respective cases.

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Directions of the product

Case Inspection: Check through a visual inspection, if there is physical damage to the case, such as scratches, dents, characteristics associated with conservation and its functionality. This technical inspection must be performed by a qualified professional.

If any part failed, it must be separated for a review and maintenance by the supplier or intended for disposal.

Packaging of instruments and / or implants in the case: The instruments and / or implants should be handled carefully and in small batches to avoid bumping or dropping the mounting assembly. All instruments and / or implants should be positioned according to the IDs of their respective cases.

After this assembly the case should be submitted to a technical inspection by an officer enabled.

Sterilization: The instruments are supplied fully decontaminated, clean and non-sterile; they must be sterilized before use. TRAUMEC recommends the sterilization methods described below.

The parameters suitable for sterilization processes (physical or chemical) for each equipment and volume must be analyzed and conducted by people trained and specialized in sterilization processes, ensuring the complete efficiency of this procedure.

In order to do so, the manufacturer instructions and methods in accordance with the internal guidelines of the hospital establishment must be followed.

In all cases, the sterilization process selected must comply with the EN556-1 standard, which establishes that the theoretical probability of vital microorganisms presence must be, at a maximum of an equal to 10⁻⁶ (S.A.L. [Sterility Assurance Level] – 10⁻⁶). The guarantee of applying a suitable sterilization process and checking if every device is sterilized are in full responsibility of the user

at every phase of the process. It is recommended that the following physical sterilization parameters in autoclaves (saturated vapor) are applied:

Cycle	Temperature	Exposure Time
Regular (pressure of 1 atm)	121°C (250°F)	30 minutes
Regular (1 atm pressure)	132°C (270°F)	15 minutes
Gravity	132°C (270°F)	45 minutes
High-Vacuum	132°C (270°F)	4 minutes

NOTE: Time must only be counted when the sterilization chamber reaches the desired temperature. The instruments must be sterilized in autoclaves before use, according to a validated method and in compliance with NBR ISO 17665 Part 1: 2010.

Reuse of Cases: The process for re-use cases for packaging implants and surgical instruments involves at least five basic steps:

Previous cleaning, decontamination, Wash, Rinse and Drying:

It is recommended that every case is cleaned immediately after surgical procedure that was employed avoiding hardening of dirt coming from the procedure.

All manual cleaning procedure must be performed using appropriate personal protective equipment. For cleaning operations on automatic equipment, the manufacturers instructions should be strictly followed, especially in caring for the product and the quality of the water used. The cases where material should be introduced open or disassembled.

The quality of water is a key factor for both the cleaning process and for the conservation of cases, therefore, the presence of particulate elements, the concentration of chemicals or substances, and pH imbalance may deteriorate the cases during the process cleaning.

Pre-cleaning: The cases must be preferably submerged and, whenever possible, opened or disassembled in a suitable container with water and a preferably enzymatic detergent at room temperature. Afterwards, it must be strictly washed under running and preferably lukewarm water. This phase must always be performed with water under 45°C, since elevated temperatures cause protein coagulation, turning the process of removing incrustations from the instrument a bit harder.

Decontamination: It is performed through submerging, whenever possible, an open or disassembled cases, in a suitable container, with a solution of disinfectant in water at room temperature (chemical disinfection), or in a warmed bath (thermochemical disinfection). The amount of time in which the instrument must be submerged depends on the operation temperature, dilution and the type of disinfectant used.

Washing: The pieces must be thoroughly brushed, with a soft bristle brush, paying special attention to articulations, grooves and gear racks. The cases, whenever relevant, must be disassembled and each component must be washed separately. Pay special attention to areas that are difficult to reach, where organic tissue, secretions or disinfection solutions may be retained or deposited.

Rinsing: The case must be plentifully rinsed under running water; articulated instruments must be opened and closed several times during the rinsing process. It is recommended to use warm water in order to rinse the case.

Drying: You must ensure that the drying processes do not introduce humidity, particles or fuzz over case surface. Pay special attention to articulations, grooves and gear racks. It is recommended that the tissue is absorbent, soft and that each part of a disassembling instrument is dried separately; and if the instrument presents cavities or slots, their inside must be completely dried.

Disposal: The disposal of disqualified parts must be performed under technical assessment and guidance. After the replacement, destroy all damaged components, preventing further improper use.

When there is the need to dispose of the instrument, it must be immediately disabled, in order to prevent its inadvertent use.

The disposal of cases must comply with the standards regarding disposal of contaminant medical waste.

We recommend that the parts are filed, bent or cut in order to be disabled. For the disposal of case, follow your local legislation in force regarding the discard of potentially contaminant products.

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6. Storage conditions

The cases must be stored in a clean, dry and ventilated place, at room temperature and away from direct light. The special storage, handling and preservation conditions of the product must be observed, in order to guarantee that the components remain intact until used in the surgical procedure. Safety measures regarding receipt, storage, transport, cleaning and batch references preservation must be adopted along the best practices of storage and distribution of medical products.

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7. Handling and transportation conditions

The case must be transported and handled as to prevent any damages or alterations of their characteristics.

The case must be managed and handled carefully, in small batches, preventing bumping or dropping. Any instruments that has been dropped or handled in an unsuitable manner or even if it is suspected that it may have suffered any damages must be isolated and submitted to the responsible technician, qualified for performing the inspection, even if it has already been through this step.

Note: Any products that has been dropped or unsuitably managed, or even if it is suspected it may have been suffered any damages ,must be sorted and segregated.

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9. Warnings

This product can be used several times since that is not damaged by wear or shock caused by handling. Therefore, every use should be made a thorough inspection, checking all aspects of the case that may affect the requirements for effectiveness and safety of the product, and when detected any irregularities discard it.

Under no circumstances, we need to use steel wool or abrasives even scouring powder to remove remaining dirt from any stage of the cleaning process.

It must be assured that every case to be sterilized is effectively cleaned. The products used must be registered at ANVISA.

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10. Precautions

The instruments must be kept in their original packages until sterilized and used. Clean it correctly after every use, in order to prevent incrustations or corrosions.

At every use, check whether the instruments have suffered any damages.

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