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**STERILE PRODUCT - PRODUCT RECOMMENDED FOR USE  
ONLY**

Application, special storage conditions, storage, handling, reuse, cautions and warnings.

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**Technical name:** Surgical Burs

**Product name:** Surgical Drills Diamond TRAUMEC

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## 1. Application

The Surgical Drills Diamond Traumec instruments are made to assist surgeons in various surgical specialties, where you have to drill, peeling, grinding bones and other tissues. These instruments are designed to be used in engines Surgical, which give full support to the use of Surgical Drills Diamond Traumec.

The Surgical Drills Diamond Traumec can be used in various applications: Mandibular osteotomy Anterior, Posterior Mandibular Osteotomy, Osteotomy Posterior Mandibular Dentoalveolar, Facial Asymmetry, Bone Reconstruction, Orbital hypertelorism, genioplasty, LeFort I osteotomy, LeFort II and LeFort III, Cantal Medial and Correction lateral Disostese mandibulofacial, Traumatic Deformities: Primary Treatment and Secondary Deformities Secondary to tumors, Apert or Crouzon syndrome and facial clefts.

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## **2. Materials**

The Surgical Drills Diamond Traumec are manufactured in AISI 420 stainless steel coated diamond in the active part.

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## **3. Form of Presentation**

The Surgical Drills Diamond Traumec are supplied sterile by ETO (Ethylene Oxide), packed in blister heat sealed, sold individually in final packaging made of triplex paper and properly labeled.

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## 5. Operating Principle

The Surgical Drills Traumec be coupled to a handpiece specifies, and this connected to a surgical motor, providing the surgeon with their modes of use: procedures for cutting bone and hard tissues. It can be used in several surgical specialties where there is the need to perform the procedures mentioned above.

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## **6. Instructions for use the product**

Opening the package by pulling the tab opening at the end of the pack, removes it from the package. Its handling should follow aseptic techniques used in operating rooms. Connect to drill in the handpiece, making sure that it is firmly attached, connect the handpiece to the surgical engine, select the power of the engine according to the needs required by the procedure and perform the operation.

The surgeon performing the procedure in charge is responsible for the choice of parameters, the selection of instruments and surgical technique required for proper execution of the surgery. As a manufacturer, Traumec does not recommend a specific surgical technique.

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## **7. Conditions of storage, shipping and handling.**

The Surgical Drills Diamond Traumec should be stored in clean, dry, airy room temperature and protected from light. The special conditions of storage, handling and storage of the product should be followed in order to ensure that the components remain intact for the procedure. Care of the receipt, storage, transportation, cleaning and maintenance of the lot of the references should be adopted in conjunction with good storage and distribution of medical products. The Surgical Drills Diamond Traumec must be transported and handled so as to prevent any damage or change in its characteristics.

The Surgical Drills Diamond Traumec must be handled carefully and individually, avoiding bumping or dropping. Any product that has been dropped or improperly handled, or suspected to have suffered damage, should be separated and segregated.

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## **8. Procedures for Re-use of Medical Product**

The process for reuse of surgical instruments involves at least five basic steps: pre-cleaning, decontamination, washing, rinsing and drying. It is recommended that all instruments be cleaned immediately after the surgical procedure that was employed to avoid the hardening of dirt taken from this procedure. Cleaning should be standardized, avoiding the spread of contamination and damage to the instrument. All cleaning procedure manual should be performed using appropriate personal protective equipment. In cleaning in automatic equipment, the manufacturers' instructions must be strictly followed, especially as regards the products and the quality of water to be used.

In no event should be used metallic brushes, steel wool or other abrasives, even the scouring, to remove dirt remaining from any stage of the cleaning process. Do not use aggressive cleaning agents, agents such as mineral acids (sulfuric, nitric). It should be ensured that the instrumental is free of any product preservation, as well as any dirt coming from the storage or repair procedure. The presence of nonsoluble products can lead to the formation of physical barriers, protecting microorganisms from the action of germicides, as well as providing the retention of unwanted dirt to further use of the instruments.

Water quality is a key factor for both the cleaning process, and for the conservation of the instruments. The presence of particulate elements, the concentration of elements or chemicals, and the imbalance of pH can damage the instrument during the cleaning process. The combination of some of these parameters may lead to fouling of precipitated minerals, can not be eliminated during removal of organic matter of scale, and the

induction process of corrosion of stainless steel, as in the case of excessive presence of chlorides. It is recommended that the water used for washing the instruments are in accordance with the quality requirements set in the sterilization process.

Note: All instruments should be cleaned immediately after the surgical procedure, thus avoids the hardening fluid originated from surgical work. The entire cleaning process should be done with extreme caution, avoiding falls, beats that could compromise the instrumental.

#### **CLEANING PRIOR:**

The Surgical Drills Diamond Traumec must be immersed in a suitable container containing water and detergent, preferably enzymatic, at room temperature. After, must be rigorously washed in water, preferably warm. This phase should be performed with water at temperatures below 45 ° C because higher temperatures cause coagulation of protein, making the removal of fouling of the instrument.

#### **DECONTAMINATION:**

The decontamination is performed by immersing the drill in a suitable container containing disinfectant solution in water at room temperature (disinfection chemistry), or in a heated bath (thermochemical disinfection). The immersion time of instrument depends on the operating temperature, the dilution, and the type of disinfectant used.

#### **WASH AND RINSE:**

The Surgical Drills Diamond Traumec should be completely brushed with soft bristle brush, paying particular attention to areas of difficult access, where it could cause the retention of tissues and deposition of secretions or disinfectant solutions and subsequently rinsed thoroughly in running water.

#### **DRYING:**

You should ensure that the drying does not cause the deposit of particles or lint on the surface of the drills. Special care should be given to areas of difficult access. It is recommended that the fabric is absorbent, soft, and each drill is dried separately, making sure that no moisture on the surface of the drills.

#### **STERILIZATION:**

In the case of reuse of Surgical Drills Diamond Traumec, we recommend the sterilization methods described below: There are adequate parameters of sterilization processes (physical or chemical) for each device and volume, which should be analyzed and conducted by trained and specialized in sterilization processes, ensuring the full effectiveness of this procedure.

To this should be followed the manufacturer's instructions and methods in accordance with internal guides in the hospital use. The sterilization process selected, in any case shall meet the standard EN556-1, which states that the theoretical probability of the presence of microorganisms vital is at most equal to 10<sup>-6</sup> (S.A.L [Sterility Assurance Level] - level of sterility assurance = 10<sup>-6</sup>). It is sole responsibility of the user to guarantee use of an appropriate sterilization process and check the sterility of all devices at any stage of the process.

It is recommended to apply the following physical parameters of sterilization in an autoclave (saturated steam):

<b>Cycle</b>	<b>Temperature</b>	<b>Exposure Time</b>
Conventional (1 atm de pressure)	121°C (250°F)	30 minutos
Conventional (1 atm de pressure)	132°C (270°F)	15 minutos
Gravity	132°C (270°F)	45 minutos
High Vacuum	132°C (270°F)	4 minutos

OBS.: The time should only be scored when the heat of the sterilization chamber reaches the desired temperature.

The Surgical Drills Diamond Traumec must be sterilized in an autoclave according to the method validated according to ISO 17665 Part 1:2010.

Disposal: The disposal of unqualified parts must be under evaluation and technical guidance. After replacing the damaged components must be destroyed in order to be avoided, inappropriately later use.

Disposal of the instrument should follow the rules relating to disposal of medical waste contaminant, dropping into bins with clear identification that it is garbage polluting.

We recommend that the pieces are cut, bent or filed for its destruction. To eliminate the instrumental follow local legal procedures of the country for disposal of potentially contaminating products.

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## 9. Warnings

Do not use if package is opened or damaged.

Do not use the drills there are signs of damage.

The perforation procedure must be carried out using the drill rotation between 500 and 1000 rpm .

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## 10. Precautions

- Never attempt to rework a drill. This will affect the specific function of its use, changing its overall dimensions. Any rework may cause excessive wear, characterizing its function, bringing great losses to the intended use and the patient.
- Continued or prolonged use of a drill can cause wear in them and result in possible disruptions. A drill worn, will inevitably longer surgical procedure, which may cause thermal necrosis due to frictional heating.
- For good surgical results and to avoid any loss of the cutting, thermal necrosis due to frictional heating, bending and risk of rupture, always use new drills every surgical procedure.
- The use of irrigation during the drilling operation will reduce the possibility of thermal necrosis.
- Ensure that the drill is properly seated and locked into the handpiece before starting engine surgery.
- Do not insert or remove the drills while the handpiece is running. Failure to follow this instruction can result in injury to staff the operating room.
- Do not apply excessive pressure, for example, bend or bow drills. Excessive pressure can bend or fracture the instruments, causing tissue damage and / or loss of tactile control. Failure to follow this instruction can result in injury to the patient and / or the operating room team.

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## 11. Contraindications.

This product can not be used in procedures in the central nervous system or central circulatory system.

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