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INSTRUCTIONS FOR USE

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TECHNICAL NAME: Box / case for Surgical Instruments

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TRADE NAME: Box Surgical instruments in Radel

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WARNING: Read carefully all of the instructions before use. Comply with all warnings and precautions mentioned in these instructions. The failure to comply with these recommendations may lead to complications.

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1. Product Identification

The boxes and trays are manufactured for packaging implants and / or instrumental used in orthopedic surgeries and bone correction, and for each type of surgery and for each set of products there is a specific case. The purpose of boxes and trays is to pack and easy transport and sterilization of instruments and materials - health products - identifying them as a way to facilitate the surgical instruments and sterilization by autoclave.

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2. Instrument Raw Materials

The following materials are employed on manufacturing of the boxes and trays: Radel

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These Instructions also regard the following products:

| CÓDIGO | DESCRIÇÃO | ILUSTRAÇÃO |
|---------------|---|---|
| PA.02.06.0006 | Tampa DIN 1/2 Radel Branco Osso Tapón DIN 1/2 Radel Blanco Hueso Cover DIN 1/2 Radel White Bone |  |
| PA.02.06.0007 | Tampa DIN 1/2 Radel Preto Tapón DIN 1/2 Radel Negro Cover DIN 1/2 Radel Black | |
| PA.02.06.0008 | Tampa DIN 1/2 Radel Amarelo Tapón DIN 1/2 Radel Amarillo Cover DIN 1/2 Radel Yellow | |
| PA.02.06.0009 | Tampa DIN 1/2 Radel Verde Tapón DIN 1/2 Radel Verde Cover DIN 1/2 Radel Green | |
| PA.02.06.0010 | Tampa DIN 1/2 Radel Azul Tapón DIN 1/2 Radel Azul Cover DIN 1/2 Radel Blue | |
| PA.02.06.0011 | Tampa DIN 1/2 Radel Vermelho Tapón DIN 1/2 Radel Rojo Cover DIN 1/2 Radel Red | |
| PA.02.06.0013 | Tampa DIN 1/2 Radel Fume Transparente Tapón DIN 1/2 Radel Gris Transparente Cover DIN 1/2 Radel Gray transparent | |
| PA.02.06.0014 | Tampa DIN 1/2 Radel Amarelo Transparente Tapón DIN 1/2 Radel Amarillo Transparente Cover DIN 1/2 Radel Yellow transparent | |
| PA.02.06.0015 | Tampa DIN 1/2 Radel Verde Transparente Tapón DIN 1/2 Radel Verde Transparente Cover DIN 1/2 Radel Green transparent | |
| PA.02.06.0016 | Tampa DIN 1/2 Radel Azul Transparente Tapón DIN 1/2 Radel Azul Transparente Cover DIN 1/2 Radel Blue Transparent | |
| PA.02.06.0017 | Tampa DIN 1/2 Radel Vermelho Transparente Tapón DIN 1/2 Radel Rojo Transparente Cover DIN 1/2 Radel Red Transparent | |

| | | |
|---------------|---|---|
| PA.02.06.0026 | Tampa DIN 3/4 Radel Branco Osso Tapón DIN 3/4 Radel Blanco Hueso Cover DIN 3/4 Radel White bone |  |
| PA.02.06.0027 | Tampa DIN 3/4 Radel Preto Tapón DIN 3/4 Radel Negro Cover DIN 3/4 Radel Black | |
| PA.02.06.0028 | Tampa DIN 3/4 Radel Amarelo Tapón DIN 3/4 Radel Amarillo Cover DIN 3/4 Radel Yellow | |
| PA.02.06.0029 | Tampa DIN 3/4 Radel Verde Tapón DIN 3/4 Radel Verde Cover DIN 3/4 Radel Green | |
| PA.02.06.0030 | Tampa DIN 3/4 Radel Azul Tapón DIN 3/4 Radel Azul Cover DIN 3/4 Radel Blue | |
| PA.02.06.0031 | Tampa DIN 3/4 Radel Vermelho Tapón DIN 3/4 Radel Rojo Cover DIN 3/4 Radel Red | |
| PA.02.06.0033 | Tampa DIN 3/4 Radel Fume Transparente Tapón DIN 3/4 Radel Gris Transparente Cover DIN 3/4 Radel Gray transparent | |
| PA.02.06.0034 | Tampa DIN 3/4 Radel Amarelo Transparente Tapón DIN 3/4 Radel Amarillo Transparente Cover DIN 3/4 Radel Yellow transparent | |
| PA.02.06.0035 | Tampa DIN 3/4 Radel Verde Transparente Tapón DIN 3/4 Radel Verde Transparente Cover DIN 3/4 Radel Green transparent | |
| PA.02.06.0036 | Tampa DIN 3/4 Radel Azul Transparente Tapón DIN 3/4 Radel Azul Transparente Cover DIN 3/4 Radel Blue transparent | |
| PA.02.06.0037 | Tampa DIN 3/4 Radel Vermelho Transparente Tapón DIN 3/4 Radel Rojo Transparente Cover DIN 3/4 Radel Red transparent | |

| | | |
|---------------|---|---|
| PA.02.06.0049 | Tampa DIN Full Radel Branco Osso Tapón DIN Full Radel Blanco Hueso Cover DIN Full Radel White Bone |  |
| PA.02.06.0050 | Tampa DIN Full Radel Preto Tapón DIN Full Radel Negro Cover DIN Full Radel Black | |
| PA.02.06.0051 | Tampa DIN Full Radel Amarelo Tapón DIN Full Radel Amarillo Cover DIN Full Radel Yellow | |
| PA.02.06.0052 | Tampa DIN Full Radel Verde Tapón DIN Full Radel Verde Cover DIN Full Radel Green | |
| PA.02.06.0053 | Tampa DIN Full Radel Azul Tapón DIN Full Radel Azul Cover DIN Full Radel Blue | |
| PA.02.06.0054 | Tampa DIN Full Radel Vermelho Tapón DIN Full Radel Rojo Cover DIN Full Radel Red | |
| PA.02.06.0056 | Tampa DIN Full Radel Fume Transparente Tapón DIN Full Radel Gris Transparente Cover DIN Full Radel Gray Transparente | |
| PA.02.06.0057 | Tampa DIN Full Radel Amarelo Transparente Tapón DIN Full Radel Amarillo Transparente Cover DIN Full Radel Yellow Transparente | |
| PA.02.06.0058 | Tampa DIN Full Radel Verde Transparente Tapón DIN Full Radel Verde Transparente Cover DIN Full Radel Green Transparente | |
| PA.02.06.0059 | Tampa DIN Full Radel Azul Transparente Tapón DIN Full Radel Azul Transparente Cover DIN Full Radel Blue Transparente | |
| PA.02.06.0060 | Tampa DIN Full Radel Vermelho Transparente Tapón DIN Full Radel Rojo Transparente Cover DIN Full Radel Red Transparente | |

| | | |
|---------------|---|---|
| PA.02.06.0061 | Base DIN 1/2 Radel Branco Osso 238x242x51 Base DIN 1/2 Radel Blanco Hueso 238x242x51 Base DIN 1/2 Radel White Bone 238x242x51 |  |
| PA.02.06.0062 | Base DIN 3/4 Radel Branco Osso 386x243x20 Base DIN 3/4 Radel Blanco Hueso 386x243x20 Base DIN 3/4 Radel White Bone 386x243x20 | |
| PA.02.06.0063 | Base DIN 3/4 Radel Branco Osso 386x243x63 Base DIN 3/4 Radel Blanco Hueso 386x243x63 Base DIN 3/4 Radel White Bone 386x243x63 | |
| PA.02.06.0064 | Base DIN Full Radel Branco Osso 514x241x38 Base DIN Full Radel Blanco Hueso 514x241x38 Base DIN Full Radel White Bone 514x241x38 | |
| PA.02.06.0065 | Base DIN Full Radel Branco Osso 514x241x76 Base DIN Full Radel Blanco Hueso 514x241x76 Base DIN Full Radel White Bone 514x241x76 | |
| PA.02.06.0066 | Base DIN Full Radel Branco Osso 514x241x101 Base DIN Full Radel Blanco Hueso 514x241x101 Base DIN Full Radel White Bone 514x241x101 | |
| PA.02.06.0067 | Bandeja 1/3 Radel Branco Osso 190x101x33 Bandeja 1/3 Radel Blanco Hueso 190x101x33 Tray 1/3 Radel White Bone 190x101x33 |  |
| PA.02.06.0068 | Bandeja 2/3 Radel Branco Osso 251x190x43 Bandeja 2/3 Radel Blanco Hueso 251x190x43 Tray 2/3 Radel White Bone 251x190x43 | |
| PA.02.06.0069 | Bandeja 2/3 Radel Branco Osso 251x190x33 Bandeja 2/3 Radel Blanco Hueso 251x190x33 Tray 2/3 Radel White Bone 251x190x33 | |
| PA.02.06.0070 | Bandeja 3/3 Radel Branco Osso 414x190x46 Bandeja 3/3 Radel Blanco Hueso 414x190x46 Tray 3/3 Radel White Bone 414x190x46 | |

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4. Operation principle

Accommodate in a systematic and organized way the implants and instruments in order to maintain the physical integrity of the instruments, facilitating handling both during surgery, as in transportation procedures, storage and sterilization.

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5. Directions

Each model of box / case for packaging surgical instruments has a specific form of packaging of instruments that are provided by the IDs each case, facilitating the recognition of the instruments. The kits do not come accompanied by instruments for surgeries, these are purchased separately, then put in their respective cases.

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PROCEDURES IN ORDER TO USE AND REUSE THE MEDICAL PRODUCT

Case Inspection: Check through a visual inspection, if there is physical damage to the case, such as scratches, dents, characteristics associated with conservation and its functionality. This technical inspection must be performed by a qualified professional.

If any part failed, it must be separated for a review and maintenance by the supplier or intended for disposal.

Packaging of instruments and / or implants in the case: The instruments and / or implants should be handled carefully and in small batches to avoid bumping or dropping the mounting assembly. All instruments and / or implants should be positioned according to the IDs of their respective cases.

After this assembly the case should be submitted to a technical inspection by an officer enabled.

Sterilization: The instruments are supplied fully decontaminated, clean and non-sterile; they must be sterilized before use. TRAUMEC recommends the sterilization methods described below.

The parameters suitable for sterilization processes (physical or chemical) for each equipment and volume must be analyzed and conducted by people trained and specialized in sterilization processes, ensuring the complete efficiency of this procedure.

In order to do so, the manufacturer instructions and methods in accordance with the internal guidelines of the hospital establishment must be followed.

In all cases, the sterilization process selected must comply with the EN556-1 standard, which establishes that the theoretical probability of vital microorganisms presence must be, at a maximum of an equal to 10⁻⁶ (S.A.L. [Sterility Assurance Level] – 10⁻⁶). The guarantee of applying a suitable sterilization process and checking if every device is sterilized are in full responsibility of the user

at every phase of the process. It is recommended that the following physical sterilization parameters in autoclaves (saturated vapor) are applied:

| Cycle | Temperature | Exposure Time |
|------------------------------|---------------|---------------|
| Regular (pressure of 1 atm) | 121°C (250°F) | 30 minutes |
| Regular (1 atm pressure) | 132°C (270°F) | 15 minutes |
| Gravity | 132°C (270°F) | 45 minutes |
| High-Vacuum | 132°C (270°F) | 4 minutes |

NOTE: Time must only be counted when the sterilization chamber reaches the desired temperature. The instruments must be sterilized in autoclaves before use, according to a validated method and in compliance with NBR ISO 17665 Part 1: 2010.

Reuse of Cases: The process for re-use cases for packaging implants and surgical instruments involves at least five basic steps:

Previous cleaning, decontamination, Wash, Rinse and Drying:

It is recommended that every case is cleaned immediately after surgical procedure that was employed avoiding hardening of dirt coming from the procedure.

All manual cleaning procedure must be performed using appropriate personal protective equipment. For cleaning operations on automatic equipment, the manufacturers instructions should be strictly followed, especially in caring for the product and the quality of the water used. The cases where material should be introduced open or disassembled.

The quality of water is a key factor for both the cleaning process and for the conservation of cases, therefore, the presence of particulate elements, the concentration of chemicals or substances, and pH imbalance may deteriorate the cases during the process cleaning.

Pre-cleaning: The cases must be preferably submerged and, whenever possible, opened or disassembled in a suitable container with water and a preferably enzymatic detergent at room temperature. Afterwards, it must be strictly washed under running and preferably lukewarm water. This phase must always be performed with water under 45°C, since elevated temperatures cause protein coagulation, turning the process of removing incrustations from the instrument a bit harder.

Decontamination: It is performed through submerging, whenever possible, an open or disassembled cases, in a suitable container, with a solution of disinfectant in water at room temperature (chemical disinfection), or in a warmed bath (thermochemical disinfection). The amount of time in which the instrument must be submerged depends on the operation temperature, dilution and the type of disinfectant used.

Washing: The pieces must be thoroughly brushed, with a soft bristle brush, paying special attention to articulations, grooves and gear racks. The cases, whenever relevant, must be disassembled and each component must be washed separately. Pay special attention to areas that are difficult to reach, where organic tissue, secretions or disinfection solutions may be retained or deposited.

Rinsing: The case must be plentifully rinsed under running water; articulated instruments must be opened and closed several times during the rinsing process. It is recommended to use warm water in order to rinse the case.

Drying: You must ensure that the drying processes do not introduce humidity, particles or fuzz over case surface. Pay special attention to articulations, grooves and gear racks. It is recommended that the tissue is absorbent, soft and that each part of a disassembling instrument is dried separately; and if the instrument presents cavities or slots, their inside must be completely dried.

Disposal: The disposal of disqualified parts must be performed under technical assessment and guidance. After the replacement, destroy all damaged components, preventing further improper use.

When there is the need to dispose of the instrument, it must be immediately disabled, in order to prevent its inadvertent use.

The disposal of cases must comply with the standards regarding disposal of contaminant medical waste.

We recommend that the parts are filed, bent or cut in order to be disabled. For the disposal of case, follow your local legislation in force regarding the discard of potentially contaminant products.

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5. Storage conditions

The box and cases must be stored in a clean, dry and ventilated place, at room temperature and away from direct light. The special storage, handling and preservation conditions of the product must be observed, in order to guarantee that the components remain intact until used in the surgical procedure. Safety measures regarding receipt, storage, transport, cleaning and batch references preservation must be adopted along the best practices of storage and distribution of medical products.

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6. Handling and transportation conditions

The case must be transported and handled as to prevent any damages or alterations of their characteristics.

The case must be managed and handled carefully, in small batches, preventing bumping or dropping. Any instruments that has been dropped or handled in an unsuitable manner or even if it is suspected that it may have suffered any damages must be isolated and submitted to the responsible technician, qualified for performing the inspection, even if it has already been through this step.

Note: Any products that has been dropped or unsuitably managed, or even if it is suspected it may have been suffered any damages ,must be sorted and segregated.

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8. Warnings

This product can be used several times since that is not damaged by wear or shock caused by handling. Therefore, every use should be made a thorough inspection, checking all aspects of the case that may affect the requirements for effectiveness and safety of the product, and when detected any irregularities discard it.

Under no circumstances, we need to use steel wool or abrasives even scouring powder to remove remaining dirt from any stage of the cleaning process.

It must be assured that every case to be sterilized is effectively cleaned. The products used must be registered at ANVISA.

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9. Precautions

The instruments must be kept in their original packages until sterilized and used. Clean it correctly after every use, in order to prevent incrustations or corrosions.

At every use, check whether the instruments have suffered any damages.

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www.traumec.com.br

Rev.: 00

Technical name: BOX / CASE FOR SURGICAL INSTRUMENTS
Trade Name: BOX SURGICAL INSTRUMENTS IN RADEL.
Business Model: COVER DIN 1/2 RADEL WHITE BONE
Code: PA.02.06.0006
Lot: XXXXXX
Quantity: 01 pc
ANVISA Registration No: 80455630024
Manufacture: XX / XX / XXXX
Validity: XX / XX / XXXX
Material: XXXXXX
Technical Director: José Luiz Caritá - CREA-SP - 0685038754

PRODUCT NOT STERILE

SPECIAL CONDITIONS OF STORAGE, STORAGE, HANDLING MEDICAL PRODUCT: SEE INSTRUCTIONS.

ANVISA Registration No 80455630024

Technical Director: José Luiz Caritá - CREA-SP - 0685038754

Manufactured by: Traumec Orthopedic Implants and Technology Import & Export Ltda.

Address: 1A JC Street, 138, Subdivision: Centennial Garden

City: Rio Claro | UF: SP | CEP -13503-510

CNPJ: 09.123.223/0001-10 - Brazilian Industry

Customer Service - Phone / Fax: +55 19 3522-1177

E-mail: sac@traumec.com.br

www.traumec.com.br

Rev.: 00